



How to Travel With Respiratory Devices

Can I bring it in my carry-on bag?



Respiratory devices such as CPAP and BiPap machines are permitted in your carry-on bag (or checked baggage).

We recommend placing the main unit in a clear plastic bag (not provided by CATSA).



Unpack the main unit from its carrying case and place it in a bin at the checkpoint.



Tubes and Masks
For sanitary reasons, they can stay in your carrying case.

▶ Attachments or parts of the machine that connect to the tubes will not be touched by the screening officer.

Exception
to the
100 ml
Rule



Distilled water for respiratory devices is considered a medical necessity and is exempt from 100 ml liquid restrictions.

Distilled water in separate containers must be presented to the screening officer for inspection.

The main unit may require **additional screening**, including explosive trace detection.

NEED MORE INFO? 1-800-OCanada (1-800-622-6232)

Visit us online! www.catsa.gc.ca/ or download our app : www.catsa.gc.ca/mobile-app